Gebr. Bräsecker GmbH & Co. KG
Troghäuser Weg 25
32657 Lemgo - Germany
Telefon +49 (0) 5261 701-400
Telefax +49 (0) 5261 701-580
info@kometmedical.de
www.kometmedical.de

Gebr. Bräsecker GmbH & Co. KG
313761 - Rev 2026-07

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Store in a dry place, ref. no. 5.3.4: Indicates that the medical device will be to be protected from moisture.

Protect from sunlight, ref. no. 5.3.2: Indicates that the medical devices will be to be protected from sources of light.

Contains hazardous substances, ref. no. 5.4.10: Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reproductive (CMR), or substances with endocrine disrupting properties. 1

Medical device, ref. no. 5.7.7: Indicates the item is a medical device.

Unique device identification, reference number 5.7.10: Indicates a carrier that contains information about a unique device identification (UDI).

Swiss Authorised Representative, reference number 5.1.2: Indicates the authorised representative in Switzerland

Regulated by European law: CE marking of conformity

Regulated by UK law: UKCA marking of conformity ("UK Conformity Assessed"), followed by the conformity assessment body's identification number

Regulated by the US Food and Drug Administration (FDA): Attention: In accordance with U.S. Federal law, this device may only be sold by or on the order of a physician (U.S. only)

Regulated by HIBC guidelines & recommendations: Healthcare Barcode

Non-regulated symbol: Quantity

Non-regulated symbol: Material used

Non-regulated symbol: GS 1 format

Material: Instruments made of corrosion resistant steel

Intended purpose:
These single-use endoprosthetic saw blades, when used in conjunction with a surgical drive system and, in some cases, a saw guide, are intended for the cutting and sectioning of bone and endogenous tissue of comparable hardness during endoprosthetic procedures. They are intended for use by surgeons and healthcare professionals in hospitals or surgical practices within an operating theatre setting on human patients, provided that the potential benefits outweighs the risk. The products themselves provide an indirect clinical benefit. The products enable work on/separation of bones as part of an treatment of the underlying condition.

Indication:
Endoprosthetics, for work on/separation of bones

Contraindications:
Stemotomy, neurosurgery.

Preventive measures/safety instructions
This medical device is marked as a single-use product and is not intended for reuse. Recommendations for safe reprocessing are therefore not available for this product. Follow the instructions on the label. Before use, check that the instrument is securely fitted. Insert the instrument into the handpiece as far as it will go. Test if the saw blade is securely seated in the drive unit while the unit is idling. Instruments/fragments can be located using X-ray imaging techniques.

Warnings:
Do not use instruments showing any kind of damage. Do not use instruments that are not adjusted so as to prevent excessive heat generation. During use (with a saw guide), avoid jamming, forcing or bending of the saw blade. During use, contact pressures must be adjusted so as to prevent damage to the instruments, malfunction, excessive vibration or excessive heat generation. Not authorized for work on metal. Do not allow the instrument to come in contact with metal during use. Risk of injury from pointed edges and sharp blades. Potential sensitization due to alloy components (nickel).

Combinations:
The operating instructions for the drive system must be followed. Only use saw blades in the saw guide if the saw blade can move freely and unimpeded. Set the saw blades in motion in the saw guide before making contact with the bone.

Storage and transport:
- transport with care
- do not use if the packaging is damaged
- protect instruments in their packaging from UV rays and high temperatures. Store in a clean and dry place. Do not store in the immediate vicinity of solvents or chemicals.

Disposal:
Dispose of instruments in tight, fracture and puncture proof containers (to prevent contamination) and, if applicable, comply with the relevant laws in effect in the country of use as well as the recommendations of the local authorities.

Procedure in case of serious adverse events
All serious adverse events occurred in connection with the device (in accordance with article 2, 58 of Regulation (EU) 2017/745) have to be reported to the manufacturer and the competent national authority.

Liability:
Applies outside the USA: The user is responsible for checking the product prior to use to ensure that it is suitable for the intended purpose. The user is responsible for the use of the instruments. In case of contributory negligence by the user, Gebr. Bräsecker partially or totally declines liability for all resulting damages, particularly if it is due to non-observance of our recommendations for use or warnings as well as the instructions for use. Applies in the USA: ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY DISCLAIMED. ALL INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOST PROFITS, INJURY TO OTHER PROPERTY, LOSS OF USE OR OTHER COMMERCIAL LOSSES, ARE EXPRESSLY DISCLAIMED. If, notwithstanding the foregoing, consequential and incidental damages cannot be excluded due to operation of law, such damages are expressly limited in amount to the purchase price of any defective products.

Brand reference:
The trademark names of third parties mentioned herein are trademarks of the respective trademark owners.

	Regulated symbols according to ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
	Catalogue number, ref. no. 5.1.6: Indicates the manufacturer's article number so that the medical device can be identified.
	Batch code, ref. no. 5.1.5: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Sterilized using irradiation, ref. no. 5.2.4: Indicates a medical device that has been sterilized using irradiation.
	Do not re-use, ref. no. 5.4.2: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Do not use if package is damaged, ref. no. 5.2.8: Indicates a medical device that should not be used if the package has been damaged or opened.
	Manufacturer, ref. no. 5.1.1: Indicates the medical device manufacturer.
	Date of manufacture, ref. no. 5.1.3: Indicates the date when the medical device was manufactured.
	Distributor, ref. no. 5.1.9: Indicates the company distributing the medical device in the region.
	Use-by date, ref. no. 5.1.4: Indicates the date after which the medical device is not to be used.
	Consult instructions for use, ref. no. 5.4.3: Indicates the need for the user to consult the instructions for use.
	Single sterile barrier system, ref. no. 5.2.11: Indicates a single sterile barrier system.
	Double sterile barrier system, ref. no. 5.2.12: Indicates two sterile barrier systems.
	Single sterile barrier system with protective packaging inside, ref. no. 5.2.13: Indicates a single sterile barrier system with protective packaging inside.
	Single sterile barrier system with protective packaging outside, ref. no. 5.2.14: Indicates a single sterile barrier system with protective packaging outside.
	Caution, ref. no. 5.4.4: Indicates the need for the user to consult the instructions for use for important safety-related information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Store in a dry place, ref. no. 5.3.4: Indicates that the medical device has to be protected from moisture.
	Protect from sunlight, ref. no. 5.3.2: Indicates that the medical devices have to be protected from sources of light.
	Contains hazardous substances, ref. no. 5.4.10: Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reproductive (CMR), or substances with endocrine disrupting properties. 1
	1 Evidence has shown that even those of the devices that contain more than 0.1 % weight percent of cobalt (CAS: 7440-48-4) do not pose a risk of carcinogenicity, mutagenicity or toxicity for reproduction (CMR classification: 1B). Therefore, no further measures are required.
	Medical device, ref. no. 5.7.7: Indicates the item is a medical device.
	Unique device identification, reference number 5.7.10: Indicates a carrier that contains information about a unique device identification (UDI).
	Swiss Authorised Representative, reference number 5.1.2: Indicates the authorised representative in Switzerland.
	Other symbols:
	Regulated by European law: CE marking of conformity
	Regulated by UK law: UKCA marking of conformity ("UK Conformity Assessed"), followed by the conformity assessment body's identification number
	Regulated by the US Food and Drug Administration (FDA): Attention: In accordance with U.S. federal law, this device may only be sold by or on the order of a physician (U.S. only)
	Regulated by HIBC guidelines & recommendations: Healthcare Barcode
	Non-regulated symbol: Quantity
	Non-regulated symbol: GS 1 format
	Material: Instruments made of corrosion resistant steel
	Intended purpose: These single-use endoprosthetic saw blades, when used in conjunction with a surgical drive system and, in some cases, a saw guide, are intended for the cutting and sectioning of bone and endogenous tissue of comparable hardness during endoprosthetic procedures. They are used exclusively by surgeons and healthcare professionals in hospitals or surgical practices within an operating theatre setting on human patients, provided that the potential benefits outweighs the risk. The products themselves provide an indirect clinical benefit. The products enable work on/separation of bones as part of the treatment of the underlying condition.
	Indication: Endoprosthetics, for work on/separation of bones
	Contraindications: Stemotomy, neurosurgery.
	Preventive measures/safety instructions This medical device is marked as a single-use product and is not intended for reuse. Recommendations for safe reprocessing are therefore not available for this product. - Follow the instructions on the label. - Before use, check that the instrument is securely fitted. - Insert the instrument into the handlepiece as far as it will go. - Test if the saw blade is securely seated in the drive unit while the unit is running. - Instruments/Fragments can be located using X-ray imaging techniques.
	Warnings: - Do not reuse instruments showing any kind of damage. - Cooling and rinsing must be adjusted so as to prevent excessive heat generation. - During use (with a saw guide), avoid jamming, leverage or bending of the saw blade. - During use, contact pressures must be adjusted so as to prevent damage to the instruments/material, excessive vibration or excessive heat generation. - Not authorized for work on metal. Do not allow the instrument to come in contact with metal during use. - Risk of injury from pointed edges and sharp blades. - Potential sensitisation due to alloy components (nickel).
	Combinations: - The operating instructions for the drive system must be followed. - Only use saw blades in the saw guide if the saw blade can move freely and unimpeded. - Set the saw blades in motion in the saw guide before making contact with the bone.
	Storage and transport: - Transport with care. - do not use if the packaging is damaged - protect instruments in their packaging from UV rays and high temperatures. Store in a clean and dry place. Do not store in the immediate vicinity of solvents or chemicals.
	Dispose: Dispose of instruments in tight, fracture and puncture proof containers (to prevent contamination) and, if applicable, comply with the relevant laws in effect in the country of use as well as the recommendations of the local authorities.
	Procedure in case of serious adverse events All serious adverse events occurred in connection with the device (in accordance with article 2, 58 of Regulation (EU) 2017 /745) have to be reported to the manufacturer and the competent national authority.
	Liability: Applies outside the USA: The user is responsible for checking the product prior to use to ensure that it is suitable for the intended purpose. The user is responsible for the use of the instrument. In case of contributory negligence by the user, Gehr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user. Applies in the USA: ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY DISCLAIMED. ALL INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOSS PROFITS, INJURY TO OTHER PROPERTY, LOSS OF USE OR OTHER COMMERCIAL LOSSES, ARE EXPRESSLY DISCLAIMED. If, notwithstanding the foregoing, consequential and incidental damages cannot be excluded due to operation of law, such damages are expressly limited in amount to the purchase price of any defective products.
	Brand reference: The trademark names of third parties mentioned herein are trademarks of the respective trademark owners.

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